## **DNaseI (CELL culture grade)**

BIONOVAS®

Catalog number :	BM0080-2000
Size :	2,000u (kunitz,25°C), 1ml
Concentration:	2u/ul
RNase activity:	None-detected
Store at -20°C	

## **Description:**

FACING TH

RESPONSIBLY

DNase I is purified from bovine pancreas that degrades single-stranded or double-stranded DNA to produce 3-hydroxyl oligonucleotides. This enzyme is used in molecular biology techniques like digestion of DNA, in the RNA isolation, nick translation and DNase I footprinting.

 Storage buffer:
 10mM HEPES (pH7.5)

 10mM CaCl2
 10mm MgCl2

 50% glycerol
 50% glycerol

 10X Reaction Buffer:
 400mM Tris-HCl (pH8.0)

 100mM MgSO4
 10mM CaCl2

Heat Inactivation (Not for RT): Add EDTA solution (pH8.0) to final 2mM, heat at 75°C for 10 min.

## **DNase I Treatment of RNA for RT-PCR**

1. Add the following components to a sterile, RNase-free tube. Keep the tube on ice during pipetting.

RNA in water 1-8µl

10X RNase-free DNase I buffer  $1\mu$ I (10X composition: 400mM Tris-HCI (pH 8.0), 100mM MgSO<sub>4</sub> and

10mM CaCl<sub>2</sub>.)

RNase-free DNase I 1u/µg RNA

Add RNase-free water to 10µl

- 2. Incubate at 37°C for 30 min.
- 3. Incubate at 70°C for 5 min to inactivate the DNase.
- 4. Add the treated RNA to the RT-PCR reaction.

## **BIONOVAS Biotechnology Co., Ltd.**

647-808-8236 No. 194 kingsdale Ave, Toronto. Ontario. M2N3W9