5X sample buffer

Product: CATALOG: Size:

FACING THE

RESPONSIBLY

5X sample buffer FA0020 1ml,25ml BIONOVAS®

Components:

0.3M Tris-HCl, pH6.8 0.25M DTT 50% Glycerol 10% SDS 0.5% NaN3 Store at -20°C(long term), 4°C (short term) Suggest aliquot and store at -20°C

BIONOVAS Biotechnology Co., Ltd.

647-808-8236 No. 194 kingsdale Ave, Toronto. Ontario. M2N3W9